

# Running workflow "trinity-transcriptom-assembly-datasets-collection"

Step 1: Input dataset collection

forward reads datasets collection

type to filter

Step 2: Input dataset collection

reverse reads datasets collection

type to filter

Step 3: FASTQ Groomer(version 1.0.4)

File to groom

Output dataset 'output' from step 1

Input FASTQ quality scores type

Sanger & Illumina 1.8+

Advanced Options

Hide Advanced Options

Step 4: FASTQ Groomer(version 1.0.4)

File to groom

Output dataset 'output' from step 2

Input FASTQ quality scores type

Sanger & Illumina 1.8+

Advanced Options

Hide Advanced Options

Step 5: Trinity(version 2.4.0.0)

Paired or Single-end data?

Paired

Left/Forward strand reads

Output dataset 'output\_file' from step 3

Right/Reverse strand reads

Output dataset 'output\_file' from step 4

Strand specific data

Not available.

Jaccard Clip options

Not available.

Run in silico normalization of reads

True

Additional Options:

Minimum Contig Length

200

Use the genome guided mode?

No

Error-corrected or circular consensus (CCS) pac bio reads

Selection is Optional ▾

Minimum count for K-mers to be assembled

1

☐ Send results to a new history

**tools used in this workflow:**

```
"tool_id": null,  
"tool_id": null,  
"tool_id": "toolshed.g2.bx.psu.edu/repos/devteam/fastq_groomer/fastq_groomer/1.0.4",  
"tool_id": "toolshed.g2.bx.psu.edu/repos/devteam/fastq_groomer/fastq_groomer/1.0.4",  
"tool_id": "toolshed.g2.bx.psu.edu/repos/iuc/trinity/trinity/2.4.0.0",
```

**tools added to the command line section in the trinity.xml file**

```
7 <command detect_errors="aggressive"><![CDATA[  
8   export PATH=/opt/galaxy/tool_dependencies/_conda/pkgs/trinity-2.4.0-5/bin:$PATH &&  
9   export PATH=/opt/galaxy/tool_dependencies/_conda/pkgs/bowtie2-2.2.6-py27_0/bin:$PATH &&  
10  export PATH=/opt/galaxy/tool_dependencies/_conda/pkgs/jellyfish-2.2.6-0/bin:$PATH &&  
11  export PATH=/opt/galaxy/tool_dependencies/_conda/pkgs/parafly-r2013_01_21-0/bin:$PATH &&
```