

Running workflow "quality-control-paired-dataset-collections"

Step 1: Input dataset collection

reads in paired dataset collection

Step 2: FASTQ Groomer(version 1.0.4)

File to groom

Output dataset 'output' from step 1

Input FASTQ quality scores type

Sanger & Illumina 1.8+

Advanced Options

Hide Advanced Options

Step 3: FastQC(version 0.65)

Short read data from your current history

Output dataset 'output_file' from step 2

Contaminant list

Selection is Optional ▾

Submodule and Limit specifying file

Selection is Optional ▾

Step 4: Trimmomatic(version 0.32.3)

Paired end data?

True

Input Type

Dataset collection pair

Select FASTQ dataset collection with R1/R2 pair

Output dataset 'output_file' from step 2

Perform initial ILLUMINACLIP step?

Not available.

Trimmomatic Operations

Trimmomatic Operation 1

Select Trimmomatic operation to perform

Sliding window trimming (SLIDINGWINDOW)

Number of bases to average across

4

Average quality required

20

Step 5: FastQC(version 0.65)

Short read data from your current history

Output dataset 'fastq_out_paired' from step 4

Contaminant list

Selection is Optional ⬆

Submodule and Limit specifying file

Selection is Optional ⬆

☐ Send results to a new history

tools used in this workflow:

```
"tool_id": null,  
"tool_id": "toolshed.g2.bx.psu.edu/repos/devteam/fastq_groomer/fastq_groomer/1.0.4",  
"tool_id": "toolshed.g2.bx.psu.edu/repos/devteam/fastqc/fastqc/0.65",  
"tool_id": "toolshed.g2.bx.psu.edu/repos/pjbriggs/trimmomatic/trimmomatic/0.32.3",  
"tool_id": "toolshed.g2.bx.psu.edu/repos/devteam/fastqc/fastqc/0.65",
```