

## Galaxy Workflow ' gatk-variant-call-pe-dataset-collection'

---

Step	Annotation
Step 1: Input dataset collection	
<b>input</b> <i>select at runtime</i>	
Step 2: Input dataset	
<b>input</b> <i>select at runtime</i>	
Step 3: FASTQ Groomer	
<b>File to groom</b> Output dataset 'output' from step 1 <b>Input FASTQ quality scores type</b> Sanger & Illumina 1.8+ <b>Advanced Options</b> Hide Advanced Options	
Step 4: Bowtie2	
<b>Is this single or paired library</b> Paired-end Dataset Collection <b>FASTQ Paired Dataset</b> Output dataset 'output_file' from step 3 <b>Write unaligned reads (in fastq format) to separate file(s)</b> False <b>Write aligned reads (in fastq format) to separate file(s)</b> False <b>Do you want to set paired-end options?</b> No <b>Will you select a reference genome from your history or use a built-in index?</b> Use a genome from the history and build index <b>Select reference genome</b> Output dataset 'output' from step 2 <b>Set read groups information?</b> Set read groups (SAM/BAM specification)	

**Auto-assign**

False

**Read group identifier (ID)**

gid

**Auto-assign**

False

**Read group sample name (SM)**

gsn

**Platform/technology used to produce the reads (PL)**

ILLUMINA

**Auto-assign**

False

**Library name (LB)**

Empty.

**Sequencing center that produced the read (CN)**

Empty.

**Description (DS)**

Empty.

**Date that run was produced (DT)**

Empty.

**Flow order (FO)**

Empty.

**The array of nucleotide bases that correspond to the key sequence of each read (KS)**

Empty.

**Programs used for processing the read group (PG)**

Empty.

**Predicted median insert size (PI)**

Not available.

**Platform unit (PU)**

Empty.

**Select analysis mode**

1: Default setting only

**Do you want to use presets?**

No, just use defaults

**Save the bowtie2 mapping statistics to the history**

False

**Step 5: Realigner Target Creator****Choose the source for the reference list**

History

**BAM file**

Output dataset 'output' from step 4

**Using reference file**

Output dataset 'output' from step 2

**Known Variants**

**Basic or Advanced GATK options**

Basic

**Basic or Advanced Analysis options**

Basic

Step 6: Indel Realigner

**Choose the source for the reference list**

History

**BAM file**

Output dataset 'output' from step 4

**Using reference file**

Output dataset 'output' from step 2

**Restrict realignment to provided intervals**

*select at runtime*

**Known Variants**

**LOD threshold above which the realigner will proceed to realign**

5.0

**Use only known indels provided as RODs**

False

**Basic or Advanced GATK options**

Basic

**Basic or Advanced Analysis options**

Basic

Step 7: Unified Genotyper

**Choose the source for the reference list**

History

**BAM files**

**BAM file 1**

**BAM file**

Output dataset 'output\_bam' from step 6

**Using reference file**

Output dataset 'output' from step 2

**Provide a dbSNP Reference-Ordered Data (ROD) file**

Don't set dbSNP

**Genotype likelihoods calculation model to employ**

BOTH

**The minimum phred-scaled confidence threshold at which variants not at 'trigger' track sites should be called**

30.0

**The minimum phred-scaled confidence threshold at which variants not at 'trigger' track sites should be emitted (and filtered if less than the calling threshold)**

30.0

**Basic or Advanced GATK options**

Basic

**Basic or Advanced Analysis options**

Basic