

# Galaxy Workflow ' quality-control-pe-dataset-collection'

Step	Annotation
------	------------

Step 1: Input dataset collection
<b>input</b> <i>select at runtime</i>

Step 2: FASTQ Groomer
<b>File to groom</b> Output dataset 'output' from step 1
<b>Input FASTQ quality scores type</b> Sanger & Illumina 1.8+
<b>Advanced Options</b> Hide Advanced Options

Step 3: FastQC
<b>Short read data from your current history</b> Output dataset 'output_file' from step 2
<b>Contaminant list</b> <i>select at runtime</i>
<b>Submodule and Limit specifying file</b> <i>select at runtime</i>

Step 4: Trimmomatic
<b>Single-end or paired-end reads?</b> Paired-end (as collection)
<b>Select FASTQ dataset collection with R1/R2 pair</b> Output dataset 'output_file' from step 2
<b>Perform initial ILLUMINACLIP step?</b> False
<b>Trimmomatic Operations</b>
<b>Trimmomatic Operation 1</b>
<b>Select Trimmomatic operation to perform</b> Sliding window trimming (SLIDINGWINDOW)
<b>Number of bases to average across</b> 4
<b>Average quality required</b> 20

Step 5: FastQC
<b>Short read data from your current history</b> Output dataset 'fastq_out_paired' from step 4

**Contaminant list***select at runtime***Submodule and Limit specifying file***select at runtime*